# COMPOUND MICROSCOPE BIMC-12A

-  [OPERATION MANUAL](#)
-  [GUIDA UTENTE](#)
-  [MANUAL DE INSTRUCCIONES](#)
-  [MANUEL D'INSTRUCTIONS](#)
-  [BEDIENUNGSANLEITUNG](#)
-  [取扱説明書](#)
-  [產品說明書](#)

1.0 Copyright Page 1 - 3

---

2.0 Microscope at a Glance Page 4

---

3.0 Introduction Page 5

---

4.0 Unpacking and Assembly Page 5

---

5.0 Microscope Accessories Page 6 - 7

---

6.0 Using the Microscope Page 8

---

7.0 Maintenance Page 9

---

8.0 Electrics Page 10

---

9.0 Recycling and Revcovery Page 10

---

## Copyright

© 2012 by JC LAB Corporation. All rights reserved. Under the copyright law, this manual cannot be reproduced in any form without the prior written permission of JC LAB. No patent liability is assumed, with respect to the use of the information contained herein.

First edition Feb. 2012

---

## Disclaimer

This manual has been validated and reviewed of accuracy. The instructions and descriptions it contains are accurate for your wifi-microscope at the time of this manual's production. However, succeeding wifi-microscope and manuals are subject to change without notice. JC LAB assumes no liability for damages incurred directly or indirectly from error, omissions or discrepancies between the wifi-microscope and the manual.

---

## Trademarks

- ◆ JC CAM logo are registered trademarks of JC LAB Corporation.
- ◆ Wi-Fi is registered trademark of Wi-Fi Alliance.
- ◆ USB, the USB logo and USB interface are trademarks or registered trademarks of USB Licensing LLC.

Other trademarks are registered trademarks not listed above may be used in this manual.

---

## EU Declaration of Conformity

This Product is carrying the CE-Mark in accordance with the related European Directives. Responsible for CE-Marking is JC LAB Japan GmbH, 5 Chome-12-13 Serigaya Kōnan-ku, Yokohama-shi, Kanagawaken 233-000, Japan. The complete and official EU Declaration of Conformity can be found on JC LAB website <http://www.jc-lab.jp> on the internet.

---

### **Disposing of the Microscope and It's Batteries**

- ◆ Discard this Microscope in accordance with applicable laws and regulations. For further information, contact your government.
  - ◆ Some models contains rechargeable batteries. After repeated use, the batteries will finally lose their ability to hold a charge and you will need to replace them. Under certain applicable laws and regulation, it may be illegal to dispose of old batteries by replacing them in the trash.
  - ◆ Please be kind to our shared environment. Check with your local government authority for details regarding where to recycle old batteries or how to dispose of them properly.
- 

### **Information for Wireless Devices**

#### **Wireless Devices and your Health**

Wireless Products, like other radio devices, emit radio frequency electromagnetic energy. The level of energy emitted by Wireless products however is far much less than the electromagnetic energy emitted by wire device like for example mobile phones.

Because Wireless products operate within the guidelines found in radio frequency safety standards and recommendations, JC LAB believes Wireless products is safe for use by consumers. These standards and recommendations reflect the consensus of scientific community and result from deliberations of panels and committees of scientists who continually review and interpret the extensive research literature.

In some situations or environments, the use of Wireless Products may be restricted by the proprietor of the building or responsible representatives of the organization. These situations may for example includes:

- ◆ Using the Wireless products equipment on board of airplanes, or
- ◆ In any other environment where the risk of interference to other devices or services is perceived or identified as harmful.

If you are uncertain of the policy that applies on the use of wireless devices in a specific organization or environment (e.g. airports), you are encouraged to ask for authorization to use the Wireless device prior turning on the equipment.

---

### **Canada — Industry Canada (IC)**

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation on this device.

The Term “IC” before the equipment certification number only signifies that the Industry Canada technical specifications were met.

---

### **USA-Federal Communications Commission (FCC)**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

---

### **JC LAB Support**

If you require any additional help your devices or if you are having problems operating devices, you may need to contact JC LAB for additional technical assistance.

#### **Before you call**

Some problem you experience may be related to software, apps or the operating system so it is importance that you investigate other sources of assistance first.

Before contacting JC LAB, try the following:

Review troubleshooting sections in the documentation supplied with your software and/or peripheral devices.

If problem occurs when you are running in your computer systems, consult the software document for troubleshooting suggestions and consider calling the system company’s technical support department for assistance.

Consult the reseller or dealer from where you purchased your devices and software— they are your best resource for current information and support.

---

#### **JC LAB technical support**


If you are still unable to solve the problems and suspect that it is hardware related, visit JC LAB technical support website

<http://www.jc-lab.jp>

---

### Microscope at a Glance

BIMC-12A

- |                            |                        |
|----------------------------|------------------------|
| 1. Monocular Head          | 11. Power Plug         |
| 2. Stand                   | 12. Fuse Box           |
| 3. Nosepiece               | 13. On/Off Switch      |
| 4. Objective Lens          | 14. Viewing adjust     |
| 5. Stage Clip              | 15. Condenser Turning  |
| 6. Mechanical Stage        | 16. Y-axis movement    |
| 7. Condenser Knob          | 17. X-axis movement    |
| 8. Fine Adjustment Knob    | 18. Light Source Cover |
| 9. Coarse Adjustment Knob  |                        |
| 10. Light Intensity Adjust |                        |

**Note: Microscope features may vary based on your location, language, carrier, and model number. Accessory parts may incur additional fees. Contact your local distributor to ensure the features.**

## 3.0 Introduction

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily classroom and laboratory use.

JC LAB reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. JC LAB declines any responsibility deriving from instrument use that does not comply with this manual.

## 4.0 Unpacking and Assembly

**4.1** The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Lift the microscope from the container and put it on a stable desk.

**4.2** The objectives are located in individual protective vials. Remove the objectives from the vials and insert them into the microscope nosepiece in the order from the lowest magnification to the highest, in a clockwise direction from the rear.

# 5.0 Microscope Accessories

**BIMC-A03A**

**Live stream with multiple devices**

BIMC-A03A can send streaming images to a maximum of 10 Wi-Fi enabled devices without the need for a router. High resolution streaming Wi-Fi camera.

**Turns conventional microscope into a wireless digital microscope camera!**

**Image Sensor:**

APTINA MT9P001 CMOS (Colorful)

**Scan Mode:**

Scan line by line

**Maximum Resolution:**

2592 x 1944 (5038848 pixels)

**Sensor Optical Format:**

1/2.5" [5.70mm x 4.28mm (H x V), Diagonal 7.13mm]

**Pixel Size:**

2.2µm x 2.2µm

**Dynamic Range:**

70.1dB

**ADC:**

12-bit, 8-bit R.G.B to PC

**SNRmax:**

38.1dB

**Spectral Characteristics:**

380-650nm (during daytime mode)

**Preview Mode (USB & Wifi 5G):**

10FPS@2592x1944, 15FPS@1920x1080, 30FPS@1280x960, 1280x720, 1024x768, 800x600, 640x480

**Binning Mode:**

1 x 1, 2 x 2

**Exposure Capability:**

3.9ms-320ms, Automatic and Manual

**White Balance:**

Automatic and Manual

**Software Interface:**

Direct Show

**Capture Modes:**

Photo and video

**Support System:**

Microsoft® Windows XP/7/8/8.1/10, IOS, Android

**Computer Configuration:**

CPU: equal or more than the second generation Intel core 2.8GHz, Memory: 2GB or more, USB port: USB 2.0 / 3.0 or more, Display: suggest 17" or larger

**Packing size:**

245 x 155 x 90mm (L x W x H)

\*Standard packing: Main Body + 0.5x Eyepiece + 30mm & 23.2mm Adapter + Software CD + USB

**Stage Micrometer for Microscope Calibration with NIST certification**

**BIMC-S09A**

**Lens Cleaning Tissue for slide cleaning**

**TOO-155E**

**Rechargeable battery (Portable power bank) For BIMC-12A**

**BIMC-ACPB12A**

**Digital Drying Cabinet, Cap. 60L**

**CHD-60**

**Slide Folder, hold 20pcs of slide**

**MCS-SL01A**

## 6.0 Using the Microscope

**6.1** Turn the observation head to a comfortable position for observation

**6.2** Ensure that the specimen is centred over the stage opening. Lock the specimen slide on the mechanical stage using the two slide clamp


**Caution:**

**It is recommended to use slides of 1.0 – 1.2 mm thickness in combination with cover glasses of 0.13 mm or 0.17 mm thickness.**

**6.3** Insert the plug of the cable into the power socket on the back of the microscope to charge it and turn the switch to the charge mode (on the right side of the main body).

**6.4** Rotate the coarse focusing knob to bring the slide into focus with the 4X objective. Adjust the fine focusing knob to get the image sharp and clear.

**6.5** Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

**Caution:**

**The maximum light intensity when using the 4x and 10x can damage the eyes!**

**6.6** Turn the nosepiece to choose the objective you need.

### 7.1 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (not condensing)
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

### 7.2 Do not!

- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

### 7.3 How to change the lamp

Before changing the lamp, you must switch it off, pull the plug out of the electrical socket. To avoid being burnt, wait until the lamp has cooled down. Unscrew the condensing lens above the lamp. Then remove the old lamp and replace it with a new one. Replace the condensing lens

**7.4** If you need to send the microscope to JC LAB for maintenance, please use the original packaging


**Universal Power:** 220 V, 50/60 Hz, DC 5V 2A

**Lamp:** 5W Cold Light s-LED

**Fuse:** T1A

## 9.0 Recycling and Recovery

Art.13 Dlsg 25 july 005 N°151. “According to directives 2002/95/EC, 2002/95/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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## ***Crafted for your world.***

At JC LAB microscope part, we take pride in the quality of our products and our user experience. That's why we offer technical advice and full post-sale support to keep your JC LAB products in optimal condition and ensure your microscopy results are of the highest quality. All of our microscopes are covered by a 5 year warranty on mechanical components and 2 year warranty on electrical components.

### **JC LAB Warranty Policy:**

All JC LAB products, mechanically and optically, are warrantied against defects in materials and workmanship under normal use and service for five years after purchase. Electrical components and digital cameras are warrantied for one year. Warranty does not cover bulbs, chargers, fuses, batteries, damage resulting from abuse or misuse, repairs or alterations performed by those other than authorized repair technicians, or damage occurring in transit. Warranty does not cover lenses that have become inoperable due to excessive dirtiness as a result of misuse or lack of normal, routine maintenance.

JC LAB must receive a completed warranty registration in order to honor any warranty repair. An official Return Material Authorization (RMA) number must accompany any return on the packing slip or on the exterior of the box. A description of the difficulty must be enclosed with the unit upon return, ship freight pre-paid to JC LAB. JC LAB will repair or replace at no charge and return. If failure was caused by misuse, alterations, accident, or abnormal conditions of operation, an estimate of repairs will be submitted for your approval prior to work being performed. If you have any questions concerning this product or warranty, contact your JC LAB local dealer

---

JC LAB Microscopes is committed to pursuing the full realization and continuous evolution of the Quality Management System in order to consolidate and improve the Company image, also through the commitment and professionalism of all Company staff.

---

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© 2012 da JC LAB Corporation. Tutti i diritti riservati. In base alla legge sul copyright, questo manuale non può essere riprodotto in alcuna forma senza previa autorizzazione scritta di JC LAB. Non si assume alcuna responsabilità di brevetto in relazione all'uso delle informazioni qui contenute. Prima edizione febbraio 2012